## Sample report

A sample report is a report containing information related to one or more samples.

## Content

A sample report contains the following content:

- 1. Table of contents
- 2. General overview of samples
- 3. Dicentric distribution related plots and tables
- 4. Filtering stats related plots and tables

## Creation

To open the sample report generation dialog, click "Report" in the menu bar at the top of the software window and select "Generate sample report".

## • Specify a name for the report

A name for this report. A report directory containing images and an report.html file will be created automatically using this name.

• Select samples

Include processed samples in the report by placing a checkmark beside appropriate samples in the "Select samples" list. Selection of one or more samples is required. Unprocessed samples are displayed but the checkbox is disabled.

- Specify SVM sigma value range Specify a range of SVM sigma values for which to generate DC distribution plots by selecting values in 'Min' and 'Max' dropdown boxes within the "Distribution of DCs in sample" area. Exclude DC distribution plots from the report if desired by unchecking the 'include' checkbox in the 'Distribution of DCs in sample' area.
- Select plots containing filtering stats Include specific filtering stats plots in the report by placing checkmarks beside appropriate plots in the 'Select plots' area.
- Click "OK" to begin report generation The report will open automatically in

	mple report				8	×
	contain informat		amples sele	ected below.		
Report Name:	SampleReport1	•				
Report Folder:	D:/Cytognomix	/ADCIDE/Repo	orts		Brows	e
Select samp	oles					^
V Healt	hCanada0Gy(73	1/731) - ISM: N	None			
🔲 Healt	hCanada 1Gy(15	65/1565) - ISM	1: None			
V Healt	hCanada2Gy(11	47/1147) - ISM	I: None			
Healt	hCanada3Gy(12	12/1212) - ISM	1: None			
V Healt	hCanada4Gy(90	9/909) - ISM: N	None			
TestS	ample1(8/8) - IS	M: None				
	ether to include	this type of pla	t hu uning		مرا الاستقالية	Ξ
	_	ma value in the Min		Max	A plot will be	
include	_			range	▼	
	_	Min		Max	▼	
Select plots	_	Min		Max		
<ul><li>✓ include</li><li>Select plots</li><li>✓ Number</li></ul>	-	Min 0.8	e specified	Max	▼	
<ul> <li>✓ include</li> <li>Select plots</li> <li>✓ Number</li> <li>✓ Mean C</li> </ul>	r of Objects	Min 0.8 gth/Width Ratio	e specified	Max	• International	
<ul> <li>✓ include</li> <li>Select plots</li> <li>✓ Number</li> <li>✓ Mean C</li> <li>✓ Mean C</li> </ul>	r of Objects hromosome Leng	Min 0.8 gth/Width Ratio	e specified	Max	A piot will be	
<ul> <li>✓ include</li> <li>Select plots</li> <li>✓ Number</li> <li>✓ Mean C</li> <li>✓ Mean C</li> <li>✓ Mean C</li> </ul>	r of Objects hromosome Leni andidate Densit	Min 0.8 gth/Width Ratio y the Difference	e specified	Max	the second	
include     Select plots     V Number     Mean C     Mean C     Mean C     Mean C     Mean C     Number	of Objects hromosome Lene andidate Densit	Min 0.8 gth/Width Ratio y te Difference Dbjects	e specified	Max	•	
include     Select plots     Vumber     Mean C     Mean C     Mean C     Number     Classifie	of Objects hromosome Len andidate Densit hromosome Finil of Segmented (	Min 0.8 gth/Width Ratio y te Difference Dbjects	e specified	Max		•

several seconds. It is also automatically saved as an HTML file located in the "Report Folder"/"Report Name" directory specified during report generation. Open the file "Report.html" to view the report. For more information about the location of saved reports consult the reports page.